## SONY

## SU-600MD

Stand specifically designed for use with Sony LCD flat panel medical monitors



#### Overview

The SU-600MD monitor stand is an optional accessory for use in healthcare environments. It's designed specifically for use with Sony LCD flat panel medical monitors.

The stylish monitor stand can be fixed to a horizontal surface or desk, or to a compatible medical device or system.

Monitor angle can be easily adjusted vertically over a -10 degree (downward) to +20 degree (upward) range by holding the top and bottom of the monitor and tilting.

The SU-600MD conforms to Video Electronics Standards Association (VESA) standards.

#### Specifications

General	
Storage/Transport Temperature	-20 °C to +60 °C
Storage/Transport Humidity	20% to 90%
Operating/Storage/Transport Pressure 700 hPa to 1060 hPa	

## SONY

Dimensions (W x H x D) *1	350 X 366.5 X 319mm
Description	Monitor Stand
Connector	VESA
Contents	Base unit (1) Neck unit (1) M5 x 20 Screws (4) M4 x 14 Screws (4)
Mass	Approx. 2.9kg (6 lb 6.3oz)
Supplied Accessories	Instructions for Use (1) Service Contact List (1)

Notes	
*1	The values for dimensions are approximate.

# SONY

### Gallery

